


Conductivity Probe Calibration Procedure with optek C8000 Converter



Make sure the ACF/ACS probe is connected correctly to the C8000 back panel !

SYSTEM SETTINGS		
DISPLAY		
DATE / TIME		
DATA LOGGER		
EVENT LOGGER		
LAMP VOLTAGE		
OPTICAL INPUTS		
SENSOR TF		
SENSOR TEMPERATURE		
18.07.2019	P01	13:32:58

SYSTEM SETTINGS		
TEMPERATURE DISPLAY	: ° CELSIUS	
SENSOR	TEMP (COND1)	
	TEMP (COND2)	
SENSOR ACTIVE	TEMP (PH1)	
	TEMP (PH2)	
CALIBRATION	: DEFAULT	
REJECT	P01	SAVE

SYSTEM SETTINGS		
TEMPERATURE DISPLAY	: ° CELSIUS	
SENSOR	: TEMP (COND1)	
SENSOR ACTIVE	NO	
	YES	
CALIBRATION	: DEFAULT	
REJECT	P01	SAVE

SYSTEM SETTINGS		
TEMPERATURE DISPLAY	: ° CELSIUS	
SENSOR	: TEMP (COND1)	
SENSOR ACTIVE	: YES	
CALIBRATION	: DEFAULT	
REJECT	Sample ACF	SAVE

Configure Temperature Sensor

- 1 From MAIN menu select **SYSTEM SETTINGS**
- 2 From SYSTEM SETTINGS select **SENSOR TEMPERATURE**
- 3 From SENSOR TEMPERATURE select **SENSOR**
- 4 Select **TEMP (COND1)** to get the temperature reading from the ACx probe at channel 1
- 5 Select **SENSOR ACTIVE**
- 6 Select **YES** to activate the temperature probe for measurement
- 7 Select **SAVE** to keep the changes

SYSTEM SETTINGS		
SENSOR TEMPERATURE		
SENSOR CONDUCTIVITY		
SENSOR pH		
SENSOR DIAGNOSTICS pH		
mA OUTPUTS		
PASSWORD PROTECTION		
PRODUCT CHANGE		
HOLD		
16.08.2019	Sample ACF	11:56:54

SYSTEM SETTINGS		
TEMPERATURE DISPLAY	:	° CELSIUS
SENSOR	:	TEMP (COND1)
SENSOR ACTIVE	:	YES
CALIBRATION	:	DEFAULT
REJECT	Sample ACF	SAVE

MAIN MENU		
MEASUREMENT DISPLAY		
DATA LOGGER DISPLAY		
SYSTEM INFORMATION DISPLAY		
PRODUCT CHANGE		
PRODUCT CONFIGURATION		
SYSTEM SETTINGS		
MAINTENANCE		
LANGUAGE / SPRACHE		
19.08.2019	Sample ACF	8:44:15

MAIN MENU		
HOLD MANUALLY		
ZERO POINT SETTING MANUALLY		
SYSTEM DIAGNOSTICS		
SENSOR ADAPTATION		
CALIBRATIONS		
SERVICE		
19.08.2019	Sample ACF	8:44:15

CALIBRATIONS		
CALIBRATION SCATTERED LIGHT		
CALIBRATION mA OUTPUTS		
CALIBRATION TEMP-INPUTS		
CALIBRATION CONDUCTIVITY MAN.		
CALIBRATION CONDUCTIVITY AUTO		
CALIBRATION pH		
CALIBRATION SUC		
16.08.2019	Sample ACF	12:05:41



CALIBRATION CONDUCTIVITY		
SENSOR : COND1		
SENSOR ACTIVE : YES		
CALIBRATION : DEFAULT		
TEMP COMPENSATION : YES		
COND (DEFAULT)	[mS/cm]:	0.9817
COND (USER)	T=24.1°C	: 0.9993
SENSOR CONSTANT : 0.9998		
OFFSET	[mS/cm]:	0.0000
REJECT	Sample ACF	SAVE

CALIBRATION CONDUCTIVITY		
SENSOR : COND1		
SENSOR ACTIVE : YES		
CALIBRATION : DEFAULT		
TEMP COMPENSATION : YES		
COND (DEFAULT)	[mS/cm]:	0.9821
COND (USER)	T=24.1°C	: 1.0000
SENSOR CONSTANT : 1.0010		
OFFSET	[mS/cm]:	0.0000
REJECT	Sample ACF	SAVE

Configure Conductivity Sensor

- 8 Select **SENSOR CONDUCTIVITY**
- 9 Select **SENSOR ACTIVE**
- 10 Select **YES** to activate the conductivity probe for measurement
- 11 Select **SAVE** to keep the changes

Configure and Perform Conductivity Calibration

- 12 Select **MAINTENANCE**
- 13 Select **CALIBRATIONS**
- 14 Select **CALIBRATION CONDUCTIVITY AUTO**
- 15 Place sensor into conductivity standard until lower part is covered with solution. The current **temperature** is displayed.
- 16 Select **COND (USER)** and enter the nominal value (mS/cm at the measured temperature), see table on bottle of conductivity solution
- 17 Wait until the display of the **SENSOR CONSTANT** is stable, then move to **SAVE** and press **ENTER**

Select User Conductivity Calibration Data for Measurement

MAIN MENU		
MEASUREMENT DISPLAY		
DATA LOGGER DISPLAY		
SYSTEM INFORMATION DISPLAY		
PRODUCT CHANGE		
PRODUCT CONFIGURATION		
SYSTEM SETTINGS		
MAINTENANCE		
LANGUAGE / SPRACHE		
19.08.2019	Sample ACF	8:51:54

18 Select **SYSTEM SETTINGS**

SYSTEM SETTINGS		
SENSOR TEMPERATURE		
SENSOR CONDUCTIVITY		
SENSOR pH		
SENSOR DIAGNOSTICS pH		
mA OUTPUTS		
PASSWORD PROTECTION		
PRODUCT CHANGE		
HOLD		
16.08.2019	Sample ACF	11:56:54

19 Select **SENSOR CONDUCTIVITY**

SYSTEM SETTINGS		
SENSOR	: COND1	
SENSOR ACTIVE	: YES	
CALIBRATION	: USER	
TEMP COMPENSATION	: YES	
REFERENCE TEMP	[°C] : 25.0000	
TEMP COEFFICIENT	[%] : 2.0000	
REJECT	Sample ACF	SAVE

20 Select **CALIBRATION** and **USER** and move to **SAVE** and confirm with **ENTER**

Verify Conductivity Calibration with Standard

CON-DATA		CON-TEMP	
mS/cm		°C	
4.0		23.4	
0.0	300.0	0.0	100.0
19.08.2019	Sample ACF		8:56:49

21 Return to **MAIN MENU** and select **MEASUREMENT DISPLAY**

22 Place the conductivity sensor into the conductivity standard solution, **to verify the calibration** is working properly.

23 The **conductivity sensor** is ready for sample measurements.



Germany

optek-Danulat GmbH
Emscherbruchallee 2
45356 Essen / Germany
Phone: +49 201 63409 0
E-Mail: info@optek.de



USA

optek-Danulat Inc.
N118 W18748 Bunsen Drive
Germantown WI 53022 / USA
Phone: +1 262 437 3600
Toll free call: +1 800 371 4288
Fax: +1 262 437 3699
E-Mail: info@optek.com



China

optek-Danulat Shanghai Co., Ltd.
Room 718 Building 1
No.88 Keyuan Road
Pudong Zhangjiang
Shanghai, China 201203
Phone: +86 21 2898 6326
Fax: +86 21 2898 6325
E-Mail: info@optek-danulat.com.cn



Singapore

optek-Danulat Pte. Ltd.
25 Int'l Business Park
#02-09 German Centre
Singapore 609916
Phone: +65 6562 8292
Fax: +65 6562 8293
E-Mail: info@optek.com.sg